



O2k and DatLab file: P ____ (A / B)		2017-		Operator:				
Sample type:		Cohort:		Sample code:				
Sample.Subsample number:		Unit:		Concentration:				
Medium: MiR05Cr								
Event	Mark name O2/Amp Slope	Mark name Amp Signal	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [µl]	A	B
			~200 µM O2		MiR			
AmR			10 µM	10		2		
HRP			1 U/ml		500 U/ml stock	4		
SOD	OBL		5 U/ml		5000 U/ml stock	2		
		1HP0	0		Mark on Amp signal			
1HP.1		1HP.1	0.1 µM	0.04	Mark on Amp signal	5		
1HP.2		1HP.2	0.2 µM	0.04	Mark on Amp signal	5		
ce	ROX				Permeabilized cells			
		2HP0	0		Mark on Amp signal			
2HP.1		2HP.1	0.1 µM	0.04	Mark on Amp signal	5		
1D	1D		1/2.5 mM	500		4/10		
1M.1	1M.1		0.1 mM	50		4		
2Oct	2Oct		0.5 mM	100		10		
2M2	2M2		2 mM	400		9.5		
3P	3P		5 mM	2000		5		
4G	4G		10 mM	2000		10		
5S	5S		50 mM	1000		100		
		3HP0	0		Mark on Amp signal			
3HP.1		3HP.1	0.1 µM	0.04	Mark on Amp signal	5		
6Gp	6Gp		10 mM	1000		20		
7U	7U		Δ0.5 µM	1	CCCP	Δ1		
8Rot	8Rot		0.5 µM	1		1		
9Ama	9Ama		2.5 µM	5		1		
		4HP0	0		Mark on Amp signal			
4HP.1		4HP.1	0.1 µM	0.04	Mark on Amp signal	5		