



<b>O2k and DatLab file: P___ ( A / B )</b>		2017-		Operator:				
Sample type:		Cohort:		Sample code:				
Sample.Subsample number:		Unit:		Concentration:				
Medium: MiR05Cr								
Event	Mark name O2/Amp Slope	Mark name Amp Signal	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [µl]	A	B
			~200 µM O2		<b>MiR</b>			
<b>AmR</b>			10 µM	10		2		
<b>HRP</b>			1 U/ml		500 U/ml stock	4		
<b>SOD</b>	<b>OBL</b>		5 U/ml		5000 U/ml stock	2		
		<b>1HP0</b>	0		<i>Mark on Amp signal</i>			
<b>1HP.1</b>		<b>1HP.1</b>	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		
<b>1HP.2</b>		<b>1HP.2</b>	0.2 µM	0.04	<i>Mark on Amp signal</i>	5		
<b>ce</b>	<b>R</b>							
<b>Dig</b>				8.1				
		<b>2HP0</b>	0		<i>Mark on Amp signal</i>			
<b>2HP.1</b>		<b>2HP.1</b>	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		
<b>1D</b>	<b>1D</b>		1/2.5 mM	500		4/10		
<b>1M.1</b>	<b>1M.1</b>		0.1 mM	50		4		
<b>2Oct</b>	<b>2Oct</b>		0.5 mM	100		10		
<b>2M2</b>	<b>2M2</b>		2 mM	400		9.5		
<b>3P</b>	<b>3P</b>		5 mM	2000		5		
<b>4G</b>	<b>4G</b>		10 mM	2000		10		
<b>5S</b>	<b>5S</b>		50 mM	1000		100		
		<b>3HP0</b>	0		<i>Mark on Amp signal</i>			
<b>3HP.1</b>		<b>3HP.1</b>	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		
<b>6Gp</b>	<b>6Gp</b>		10 mM	1000		20		
<b>7U</b>	<b>7U</b>		Δ0.5 µM	1	CCCP	Δ1		
<b>8Rot</b>	<b>8Rot</b>		0.5 µM	1		1		
<b>9Ama</b>	<b>9Ama</b>		2.5 µM	5		1		
		<b>4HP0</b>	0		<i>Mark on Amp signal</i>			
<b>4HP.1</b>		<b>4HP.1</b>	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		